

Antibodies Product Datasheet

Product Name:	GBA Antibody
Catalog Number:	C50794
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GBA
Concentration:	1 mg/ml
Gene Name:	GBA GC GLUC
Protein Name:	Glucosylceramidase (EC 3.2.1.45) (Acid beta-glucosidase) (Alglucerase) (Beta-glucocerebrosidase) (D-glucosyl-N-acylsphingosine glucosylhydrolase) (Imiglucerase)
Other Name:	
Human Swiss Prot No.:	P04062
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GBA at Human, Mouse
Storage and Stability:	Store at -20° C/1 year