

Antibodies Product Datasheet

Product Name:	ACSM2A Antibody
Catalog Number:	C50828
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ACSM2A
Concentration:	1 mg/ml
Gene Name:	ACSM2A ACSM2 MACS2
Protein Name:	Acyl-coenzyme A synthetase ACSM2A, mitochondrial (EC 6.2.1.2) (Acyl-CoA synthetase medium-chain family member 2A) (Butyrate--CoA ligase 2A) (Butyryl-coenzyme A synthetase 2A) (Middle-chain acyl-CoA synthetase 2A)
Other Name:	
Human Swiss Prot No.:	Q08AH3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ACSM2A at Human
Storage and Stability:	Store at -20° C/1 year