

Antibodies Product Datasheet

Product Name:	IGBP1 Antibody
Catalog Number:	C50837
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human IGBP1
Concentration:	1 mg/ml
Gene Name:	IGBP1 IBP1
Protein Name:	Immunoglobulin-binding protein 1 (B-cell signal transduction molecule alpha 4) (Protein alpha-4) (CD79a-binding protein 1) (Protein phosphatase 2/4/6 regulatory subunit) (Renal carcinoma antigen NY-REN-16)
Other Name:	
Human Swiss Prot No.:	P78318
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IGBP1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year