

Antibodies Product Datasheet

Product Name:	LRCH1 Antibody
Catalog Number:	C50845
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human LRCH1
Concentration:	1 mg/ml
Gene Name:	LRCH1 CHDC1 KIAA1016
Protein Name:	Leucine-rich repeat and calponin homology domain-containing protein 1 (Calponin homology domain-containing protein 1) (Neuronal protein 81) (NP81)
Other Name:	
Human Swiss Prot No.:	Q9Y2L9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LRCH1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year