

---

## Antibodies Product Datasheet

|                               |   |
|-------------------------------|---|
| <b>Product Name:</b>          | 1A6/DRIM Antibody   |
| <b>Catalog Number:</b>        | C50853  |
| <b>Reactivity:</b>            | Human   |
| <b>Applications:</b>          | WB  |
| <b>Dilutions:</b>             | WB 1:500-2000   |
| <b>Modification:</b>          | Unmodified/Total  |
| <b>Source:</b>                | Polyclonal, Rabbit, IgG   |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                       |
| <b>Immunogen:</b>             | Synthesized peptide derived from human 1A6/DRIM   |
| <b>Concentration:</b>         | 1 mg/ml   |
| <b>Gene Name:</b>             | UTP20 DRIM  |
| <b>Protein Name:</b>          | Small subunit processome component 20 homolog (Down-regulated in metastasis protein) (Novel nucleolar protein 73) (NNP73) (Protein Key-1A6) |
| <b>Other Name:</b>            |   |
| <b>Human Swiss Prot No.:</b>  | O75691  |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.   |
| <b>Specificity:</b>           | This antibody detects endogenous levels of 1A6/DRIM at Human  |
| <b>Storage and Stability:</b> | Store at -20° C/1 year  |