

Antibodies Product Datasheet

Product Name:	ANKS1B Antibody
Catalog Number:	C50882
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ANKS1B
Concentration:	1 mg/ml
Gene Name:	ANKS1B
Protein Name:	Ankyrin repeat and sterile alpha motif domain-containing protein 1B (Amyloid-beta protein intracellular domain-associated protein 1) (AIDA-1) (E2A-PBX1-associated protein) (EB-1)
Other Name:	
Human Swiss Prot No.:	Q7Z6G8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ANKS1B at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year