

Antibodies Product Datasheet

Product Name:	Kir2.2 Antibody
Catalog Number:	C50885
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human Kir2.2
Concentration:	1 mg/ml
Gene Name:	KCNJ12 IRK2 KCNJN1
Protein Name:	ATP-sensitive inward rectifier potassium channel 12 (Inward rectifier K(+) channel Kir2.2) (IRK-2) (Inward rectifier K(+) channel Kir2.2v) (Potassium channel, inwardly rectifying subfamily J member 12)
Other Name:	
Human Swiss Prot No.:	Q14500
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Kir2.2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year