

Antibodies Product Datasheet

Product Name:	DHDH Antibody
Catalog Number:	C50894
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DHDH
Concentration:	1 mg/ml
Gene Name:	DHDH 2DD
Protein Name:	Trans-1,2-dihydrobenzene-1,2-diol dehydrogenase (EC 1.3.1.20) (D-xylose 1-dehydrogenase) (D-xylose-NADP dehydrogenase) (EC 1.1.1.179) (Dimeric dihydrodiol dehydrogenase) (Hum2DD)
Other Name:	
Human Swiss Prot No.:	Q9UQ10
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DHDH at Human, Mouse
Storage and Stability:	Store at -20° C/1 year