

## Antibodies Product Datasheet

|                               |                                                                                                                       |
|-------------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | HS2ST1 Antibody                                                                                                       |
| <b>Catalog Number:</b>        | C50895                                                                                                                |
| <b>Reactivity:</b>            | Human, Mouse                                                                                                          |
| <b>Applications:</b>          | WB                                                                                                                    |
| <b>Dilutions:</b>             | WB 1:500-2000                                                                                                         |
| <b>Modification:</b>          | Unmodified/Total                                                                                                      |
| <b>Source:</b>                | Polyclonal, Rabbit,IgG                                                                                                |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Immunogen:</b>             | Synthesized peptide derived from human HS2ST1                                                                         |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                               |
| <b>Gene Name:</b>             | HS2ST1 HS2ST KIAA0448                                                                                                 |
| <b>Protein Name:</b>          | Heparan sulfate 2-O-sulfotransferase 1 (2-O-sulfotransferase) (2OST)<br>(EC 2.8.2.-)                                  |
| <b>Other Name:</b>            |                                                                                                                       |
| <b>Human Swiss Prot No.:</b>  | Q7LGA3                                                                                                                |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| <b>Specificity:</b>           | This antibody detects endogenous levels of HS2ST1 at Human, Mouse                                                     |
| <b>Storage and Stability:</b> | Store at -20° C/1 year                                                                                                |