

Antibodies Product Datasheet

Product Name:	ALG11 Antibody
Catalog Number:	C50907
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ALG11
Concentration:	1 mg/ml
Gene Name:	ALG11 GT8
Protein Name:	GDP-Man:Man(3)GlcNAc(2)-PP-Dol alpha-1,2-mannosyltransferase (EC 2.4.1.131) (Asparagine-linked glycosylation protein 11 homolog) (Glycolipid 2-alpha-mannosyltransferase)
Other Name:	
Human Swiss Prot No.:	Q2TAA5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ALG11 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year