

Antibodies Product Datasheet

Product Name:	DMWD Antibody
Catalog Number:	C50914
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DMWD
Concentration:	1 mg/ml
Gene Name:	DMWD DM9
Protein Name:	Dystrophia myotonica WD repeat-containing protein (Dystrophia myotonica-containing WD repeat motif protein) (Protein 59) (Protein DMR-N9)
Other Name:	
Human Swiss Prot No.:	Q09019
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DMWD at Human, Mouse
Storage and Stability:	Store at -20° C/1 year