

Antibodies Product Datasheet

Product Name:	GALNT7 Antibody
Catalog Number:	C50917
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GALNT7
Concentration:	1 mg/ml
Gene Name:	GALNT7
Protein Name:	N-acetylgalactosaminyltransferase 7 (EC 2.4.1.-) (Polypeptide GalNAc transferase 7) (GalNAc-T7) (pp-GaNTase 7) (Protein-UDP acetylgalactosaminyltransferase 7) (UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase 7)
Other Name:	
Human Swiss Prot No.:	Q86SF2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GALNT7 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year