

Antibodies Product Datasheet

Product Name:	FKBP6 Antibody
Catalog Number:	C50931
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FKBP6
Concentration:	1 mg/ml
Gene Name:	FKBP6 FKBP36
Protein Name:	Peptidyl-prolyl cis-trans isomerase FKBP6 (PPIase FKBP6) (EC 5.2.1.8) (36 kDa FK506-binding protein) (36 kDa FKBP) (FKBP-36) (FK506-binding protein 6) (FKBP-6) (Immunophilin FKBP36) (Rotamase)
Other Name:	
Human Swiss Prot No.:	O75344
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FKBP6 at Human
Storage and Stability:	Store at -20° C/1 year