

Antibodies Product Datasheet

Product Name:	GTF2F2 Antibody
Catalog Number:	C50938
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GTF2F2
Concentration:	1 mg/ml
Gene Name:	GTF2F2 RAP30
Protein Name:	General transcription factor IIF subunit 2 (EC 3.6.4.12) (ATP-dependent helicase GTF2F2) (General transcription factor IIF 30 kDa subunit) (Transcription initiation factor IIF subunit beta) (TFIIF-beta) (Transcription initiation factor RAP30)
Other Name:	
Human Swiss Prot No.:	P13984
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GTF2F2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year