

Antibodies Product Datasheet

Product Name:	BRUNOL5 Antibody
Catalog Number:	C50942
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human BRUNOL5
Concentration:	1 mg/ml
Gene Name:	CELF5 BRUNOL5
Protein Name:	CUGBP Elav-like family member 5 (CELF-5) (Bruno-like protein 5) (CUG-BP- and ETR-3-like factor 5) (RNA-binding protein BRUNOL-5)
Other Name:	
Human Swiss Prot No.:	Q8N6W0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BRUNOL5 at Human
Storage and Stability:	Store at -20° C/1 year