

Antibodies Product Datasheet

Product Name:	ETF1 Antibody
Catalog Number:	C50952
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ETF1
Concentration:	1 mg/ml
Gene Name:	ETF1 ERF1 RF1 SUP45L1
Protein Name:	Eukaryotic peptide chain release factor subunit 1 (Eukaryotic release factor 1) (eRF1) (Protein Cl1) (TB3-1)
Other Name:	
Human Swiss Prot No.:	P62495
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ETF1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year