

Antibodies Product Datasheet

Product Name:	HDDC2 Antibody
Catalog Number:	C50959
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HDDC2
Concentration:	1 mg/ml
Gene Name:	HDDC2 C6orf74 NS5ATP2 CGI-130
Protein Name:	HD domain-containing protein 2 (Hepatitis C virus NS5A-transactivated protein 2) (HCV NS5A-transactivated protein 2)
Other Name:	
Human Swiss Prot No.:	Q7Z4H3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HDDC2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year