

Antibodies Product Datasheet

Product Name:	FOX2/RBM9 Antibody
Catalog Number:	C50973
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FOX2/RBM9
Concentration:	1 mg/ml
Gene Name:	RBFOX2 FOX2 HRNBP2 RBM9 RTA
Protein Name:	RNA binding protein fox-1 homolog 2 (Fox-1 homolog B) (Hexaribonucleotide-binding protein 2) (RNA-binding motif protein 9) (RNA-binding protein 9) (Repressor of tamoxifen transcriptional activity)
Other Name:	
Human Swiss Prot No.:	O43251
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FOX2/RBM9 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year