

Antibodies Product Datasheet

Product Name:	CCT6B Antibody
Catalog Number:	C51013
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CCT6B
Concentration:	1 mg/ml
Gene Name:	CCT6B
Protein Name:	T-complex protein 1 subunit zeta-2 (TCP-1-zeta-2) (CCT-zeta-2) (CCT-zeta-like) (TCP-1-zeta-like) (Testis-specific Tcp20) (Testis-specific protein TSA303)
Other Name:	
Human Swiss Prot No.:	Q92526
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CCT6B at Human, Mouse
Storage and Stability:	Store at -20° C/1 year