

Antibodies Product Datasheet

Product Name:	DHX37 Antibody
Catalog Number:	C51019
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DHX37
Concentration:	1 mg/ml
Gene Name:	DHX37 DDX37 KIAA1517
Protein Name:	Probable ATP-dependent RNA helicase DHX37 (EC 3.6.4.13) (DEAH box protein 37)
Other Name:	
Human Swiss Prot No.:	Q8IY37
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DHX37 at Human
Storage and Stability:	Store at -20° C/1 year