

Antibodies Product Datasheet

Product Name:	KDM6B/JMJD3 Antibody
Catalog Number:	C51049
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KDM6B/JMJD3
Concentration:	1 mg/ml
Gene Name:	KDM6B JMJD3 KIAA0346
Protein Name:	Lysine-specific demethylase 6B (EC 1.14.11.-) (JmjC domain-containing protein 3) (Jumonji domain-containing protein 3) (Lysine demethylase 6B)
Other Name:	
Human Swiss Prot No.:	O15054
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KDM6B/JMJD3 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year