

Antibodies Product Datasheet

Product Name:	Hepcidin-25 Antibody
Catalog Number:	C51058
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human Hepcidin-25
Concentration:	1 mg/ml
Gene Name:	HAMP HEPC LEAP1 UNQ487/PRO1003
Protein Name:	Hepcidin (Liver-expressed antimicrobial peptide 1) (LEAP-1) (Putative liver tumor regressor) (PLTR) [Cleaved into: Hepcidin-25 (Hepc25)Hepcidin-20 (Hepc20)]
Other Name:	
Human Swiss Prot No.:	P81172
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Hepcidin-25 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year