

Antibodies Product Datasheet

Product Name:	INPP4B Antibody
Catalog Number:	C51077
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human INPP4B
Concentration:	1 mg/ml
Gene Name:	INPP4B
Protein Name:	Type II inositol 3,4-bisphosphate 4-phosphatase (EC 3.1.3.66) (Inositol polyphosphate 4-phosphatase type II)
Other Name:	
Human Swiss Prot No.:	O15327
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of INPP4B at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year