

Antibodies Product Datasheet

Product Name:	CARD15/NOD2 Antibody
Catalog Number:	C51081
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CARD15/NOD2
Concentration:	1 mg/ml
Gene Name:	NOD2 CARD15 IBD1
Protein Name:	Nucleotide-binding oligomerization domain-containing protein 2 (Caspase recruitment domain-containing protein 15) (Inflammatory bowel disease protein 1)
Other Name:	
Human Swiss Prot No.:	Q9HC29
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CARD15/NOD2 at Human
Storage and Stability:	Store at -20° C/1 year