

Antibodies Product Datasheet

Product Name:	Leptin Receptor (Phospho Tyr986) Antibody
Catalog Number:	C51085
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Leptin Receptor (Phospho Tyr986)
Concentration:	1 mg/ml
Gene Name:	LEPR DB OBR
Protein Name:	Leptin receptor (LEP-R) (HuB219) (OB receptor) (OB-R) (CD antigen CD295)
Other Name:	Leptin receptor (LEP-R) (HuB219) (OB receptor) (OB-R) (CD antigen CD295)
Human Swiss Prot No.:	P48357
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Leptin Receptor (Phospho Tyr986) Rabbit pAb at Human, Mouse, Rat
Storage and Stability:	Store at -20° C/1 year