

## Antibodies Product Datasheet

|                               |                                                                                                                                                                   |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | Bcl-6 (Phospho Ser343) Antibody                                                                                                                                   |
| <b>Catalog Number:</b>        | C51087                                                                                                                                                            |
| <b>Reactivity:</b>            | Human, Mouse                                                                                                                                                      |
| <b>Applications:</b>          | IHC, WB                                                                                                                                                           |
| <b>Dilutions:</b>             | WB 1:500-2000 IHC 1:50-200                                                                                                                                        |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                                                  |
| <b>Source:</b>                | Rabbit, polyclonal                                                                                                                                                |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.                                                         |
| <b>Immunogen:</b>             | Synthesized peptide derived from human Bcl-6 (Phospho Ser343)                                                                                                     |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                           |
| <b>Gene Name:</b>             | BCL6 BCL5 LAZ3 ZBTB27 ZNF51                                                                                                                                       |
| <b>Protein Name:</b>          | B-cell lymphoma 6 protein (BCL-6) (B-cell lymphoma 5 protein) (BCL-5) (Protein LAZ-3) (Zinc finger and BTB domain-containing protein 27) (Zinc finger protein 51) |
| <b>Other Name:</b>            | B-cell lymphoma 6 protein (BCL-6) (B-cell lymphoma 5 protein) (BCL-5) (Protein LAZ-3) (Zinc finger and BTB domain-containing protein 27) (Zinc finger protein 51) |
| <b>Human Swiss Prot No.:</b>  | P41182                                                                                                                                                            |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.                                                                                                    |
| <b>Specificity:</b>           | This antibody detects endogenous levels of Bcl-6 (Phospho Ser343) Rabbit pAb at Human, Mouse                                                                      |
| <b>Storage and Stability:</b> | Store at -20° C/1 year                                                                                                                                            |