

Antibodies Product Datasheet

Product Name:	C-RAF (Phospho Tyr340) Antibody
Catalog Number:	C51088
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human C-RAF (Phospho Tyr340)
Concentration:	1 mg/ml
Gene Name:	RAF1 RAF
Protein Name:	RAF proto-oncogene serine/threonine-protein kinase (EC 2.7.11.1) (Proto-oncogene c-RAF) (cRaf) (Raf-1)
Other Name:	RAF proto-oncogene serine/threonine-protein kinase (EC 2.7.11.1) (Proto-oncogene c-RAF) (cRaf) (Raf-1)
Human Swiss Prot No.:	P04049
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of C-RAF (Phospho Tyr340) Rabbit pAb at Human, Mouse, Rat
Storage and Stability:	Store at -20° C/1 year