

Antibodies Product Datasheet

Product Name:	Estrogen Receptor- α (Phospho Thr311) Antibody
Catalog Number:	C51090
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Estrogen Receptor- α (Phospho Thr311)
Concentration:	1 mg/ml
Gene Name:	ESR1 ESR NR3A1
Protein Name:	Estrogen receptor (ER) (ER-alpha) (Estradiol receptor) (Nuclear receptor subfamily 3 group A member 1)
Other Name:	Estrogen receptor (ER) (ER-alpha) (Estradiol receptor) (Nuclear receptor subfamily 3 group A member 1)
Human Swiss Prot No.:	P03372
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Estrogen Receptor- α (Phospho Thr311) Rabbit pAb at Human, Mouse, Rat
Storage and Stability:	Store at -20° C/1 year