

Antibodies Product Datasheet

Product Name:	CTNND1 (Phospho Tyr257) Antibody
Catalog Number:	C51093
Reactivity:	Human, Mouse
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CTNND1 (Phospho Tyr257)
Concentration:	1 mg/ml
Gene Name:	CTNND1 KIAA0384
Protein Name:	Catenin delta-1 (Cadherin-associated Src substrate) (CAS) (p120 catenin) (p120(ctn)) (p120(cas))
Other Name:	Catenin delta-1 (Cadherin-associated Src substrate) (CAS) (p120 catenin) (p120(ctn)) (p120(cas))
Human Swiss Prot No.:	O60716
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CTNND1 (Phospho Tyr257) Rabbit pAb at Human, Mouse
Storage and Stability:	Store at -20° C/1 year