

Antibodies Product Datasheet

Product Name:	MYPT1/2 (Phospho Ser668/Ser618) Antibody
Catalog Number:	C51097
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MYPT1/2 (Phospho Ser668/Ser618)
Concentration:	1 mg/ml
Gene Name:	PPP1R12A MBS MYPT1 MYPT2
Protein Name:	Protein phosphatase 1 regulatory subunit 12A (Myosin phosphatase-targeting subunit 1) (Myosin phosphatase target subunit 1) (Protein phosphatase myosin-binding subunit)
Other Name:	Protein phosphatase 1 regulatory subunit 12A (Myosin phosphatase-targeting subunit 1) (Myosin phosphatase target subunit 1) (Protein phosphatase myosin-binding subunit)
Human Swiss Prot No.:	O14974
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of MYPT1/2 (Phospho Ser668/Ser618) Rabbit pAb at Human, Mouse, Rat
Storage and Stability:	Store at -20° C/1 year