

Antibodies Product Datasheet

Product Name:	ATM (Phospho Ser367) Antibody
Catalog Number:	C51101
Reactivity:	Human, Mouse
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ATM (Phospho Ser367)
Concentration:	1 mg/ml
Gene Name:	ATM
Protein Name:	Serine-protein kinase ATM (EC 2.7.11.1) (Ataxia telangiectasia mutated) (A-T mutated)
Other Name:	Serine-protein kinase ATM (EC 2.7.11.1) (Ataxia telangiectasia mutated) (A-T mutated)
Human Swiss Prot No.:	Q13315
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ATM (Phospho Ser367) Rabbit pAb at Human, Mouse
Storage and Stability:	Store at -20° C/1 year