

Antibodies Product Datasheet

Product Name:	BRCA1 (Phospho Ser1497) Antibody
Catalog Number:	C51104
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human BRCA1 (Phospho Ser1497)
Concentration:	1 mg/ml
Gene Name:	BRCA1 RNF53
Protein Name:	Breast cancer type 1 susceptibility protein (EC 6.3.2.-) (RING finger protein 53)
Other Name:	Breast cancer type 1 susceptibility protein (EC 6.3.2.-) (RING finger protein 53)
Human Swiss Prot No.:	P38398
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BRCA1 (Phospho Ser1497) Rabbit pAb at Human, Mouse, Rat
Storage and Stability:	Store at -20° C/1 year