

Antibodies Product Datasheet

Product Name:	ALOX5 (Phospho Ser663) Antibody
Catalog Number:	C51107
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ALOX5 (Phospho Ser663)
Concentration:	1 mg/ml
Gene Name:	ALOX5 LOG5
Protein Name:	Arachidonate 5-lipoxygenase (5-LO) (5-lipoxygenase) (EC 1.13.11.34)
Other Name:	Arachidonate 5-lipoxygenase (5-LO) (5-lipoxygenase) (EC 1.13.11.34)
Human Swiss Prot No.:	P09917
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ALOX5 (Phospho Ser663) Rabbit pAb at Human, Mouse, Rat
Storage and Stability:	Store at -20° C/1 year