

Antibodies Product Datasheet

Product Name:	EPHA1 (Phospho Tyr781) Antibody
Catalog Number:	C51111
Reactivity:	Human, Mouse
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EPHA1 (Phospho Tyr781)
Concentration:	1 mg/ml
Gene Name:	EPHA1 EPH EPHT EPHT1
Protein Name:	Ephrin type-A receptor 1 (hEpha1) (EC 2.7.10.1) (EPH tyrosine kinase) (EPH tyrosine kinase 1) (Erythropoietin-producing hepatoma receptor) (Tyrosine-protein kinase receptor EPH)
Other Name:	Ephrin type-A receptor 1 (hEpha1) (EC 2.7.10.1) (EPH tyrosine kinase) (EPH tyrosine kinase 1) (Erythropoietin-producing hepatoma receptor) (Tyrosine-protein kinase receptor EPH)
Human Swiss Prot No.:	P21709
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of EPHA1 (Phospho Tyr781) Rabbit pAb at Human, Mouse
Storage and Stability:	Store at -20° C/1 year