

Antibodies Product Datasheet

Product Name:	BTK (Phospho Ser179) Antibody
Catalog Number:	C51114
Reactivity:	Human, Mouse
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human BTK (Phospho Ser179)
Concentration:	1 mg/ml
Gene Name:	BTK AGMX1 ATK BPK
Protein Name:	Tyrosine-protein kinase BTK (EC 2.7.10.2) (Agammaglobulinaemia tyrosine kinase) (ATK) (B-cell progenitor kinase) (BPK) (Bruton tyrosine kinase)
Other Name:	Tyrosine-protein kinase BTK (EC 2.7.10.2) (Agammaglobulinaemia tyrosine kinase) (ATK) (B-cell progenitor kinase) (BPK) (Bruton tyrosine kinase)
Human Swiss Prot No.:	Q06187
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BTK (Phospho Ser179) Rabbit pAb at Human, Mouse
Storage and Stability:	Store at -20° C/1 year