

Antibodies Product Datasheet

Product Name:	FANCD2 (Phospho Thr691) Antibody
Catalog Number:	C51116
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FANCD2 (Phospho Thr691)
Concentration:	1 mg/ml
Gene Name:	FANCD2 FACD
Protein Name:	Fanconi anemia group D2 protein (Protein FACD2)
Other Name:	Fanconi anemia group D2 protein (Protein FACD2)
Human Swiss Prot No.:	Q9BXW9
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FANCD2 (Phospho Thr691) Rabbit pAb at Human, Mouse, Rat
Storage and Stability:	Store at -20° C/1 year