

Antibodies Product Datasheet

Product Name:	eNOS (Phospho Ser114) Antibody
Catalog Number:	C51120
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human eNOS (Phospho Ser114)
Concentration:	1 mg/ml
Gene Name:	NOS3
Protein Name:	Nitric oxide synthase, endothelial (EC 1.14.13.39) (Constitutive NOS) (cNOS) (EC-NOS) (Endothelial NOS) (eNOS) (NOS type III) (NOSIII)
Other Name:	Nitric oxide synthase, endothelial (EC 1.14.13.39) (Constitutive NOS) (cNOS) (EC-NOS) (Endothelial NOS) (eNOS) (NOS type III) (NOSIII)
Human Swiss Prot No.:	P29474
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of eNOS (Phospho Ser114) Rabbit pAb at Human, Mouse, Rat
Storage and Stability:	Store at -20° C/1 year