

Antibodies Product Datasheet

Product Name:	FKHRL1 (Phospho Ser321) Antibody
Catalog Number:	C51124
Reactivity:	Human, Mouse
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FKHRL1 (Phospho Ser321)
Concentration:	1 mg/ml
Gene Name:	FOXO3 FKHRL1 FOXO3A
Protein Name:	Forkhead box protein O3 (AF6q21 protein) (Forkhead in rhabdomyosarcoma-like 1)
Other Name:	Forkhead box protein O3 (AF6q21 protein) (Forkhead in rhabdomyosarcoma-like 1)
Human Swiss Prot No.:	O43524
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FKHRL1 (Phospho Ser321) Rabbit pAb at Human, Mouse
Storage and Stability:	Store at -20° C/1 year