

Antibodies Product Datasheet

Product Name:	Chk2 (Phospho Ser379) Antibody
Catalog Number:	C51129
Reactivity:	Human, Mouse
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Chk2 (Phospho Ser379)
Concentration:	1 mg/ml
Gene Name:	CHEK2 CDS1 CHK2 RAD53
Protein Name:	Serine/threonine-protein kinase Chk2 (EC 2.7.11.1) (CHK2 checkpoint homolog) (Cds1 homolog) (Hucds1) (hCds1) (Checkpoint kinase 2)
Other Name:	Serine/threonine-protein kinase Chk2 (EC 2.7.11.1) (CHK2 checkpoint homolog) (Cds1 homolog) (Hucds1) (hCds1) (Checkpoint kinase 2)
Human Swiss Prot No.:	O96017
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Chk2 (Phospho Ser379) Rabbit pAb at Human, Mouse
Storage and Stability:	Store at -20° C/1 year