

Antibodies Product Datasheet

Product Name:	Connexin 43 (Phospho Ser282) Antibody
Catalog Number:	C51130
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Connexin 43 (Phospho Ser282)
Concentration:	1 mg/ml
Gene Name:	GJA1 GJAL
Protein Name:	Gap junction alpha-1 protein (Connexin-43) (Cx43) (Gap junction 43 kDa heart protein)
Other Name:	Gap junction alpha-1 protein (Connexin-43) (Cx43) (Gap junction 43 kDa heart protein)
Human Swiss Prot No.:	P17302
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Connexin 43 (Phospho Ser282) Rabbit pAb at Human, Mouse, Rat
Storage and Stability:	Store at -20° C/1 year