

Antibodies Product Datasheet

Product Name:	CSK (Phospho Tyr304) Antibody
Catalog Number:	C51134
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CSK (Phospho Tyr304)
Concentration:	1 mg/ml
Gene Name:	CSK
Protein Name:	Tyrosine-protein kinase CSK (EC 2.7.10.2) (C-Src kinase) (Protein-tyrosine kinase CYL)
Other Name:	Tyrosine-protein kinase CSK (EC 2.7.10.2) (C-Src kinase) (Protein-tyrosine kinase CYL)
Human Swiss Prot No.:	P41240
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CSK (Phospho Tyr304) Rabbit pAb at Human, Mouse, Rat
Storage and Stability:	Store at -20° C/1 year