

Antibodies Product Datasheet

Product Name:	Bim (Phospho Thr56/116) Antibody
Catalog Number:	C51136
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Bim (Phospho Thr56/116)
Concentration:	1 mg/ml
Gene Name:	BCL2L11 BIM
Protein Name:	Bcl-2-like protein 11 (Bcl2-L-11) (Bcl2-interacting mediator of cell death)
Other Name:	Bcl-2-like protein 11 (Bcl2-L-11) (Bcl2-interacting mediator of cell death)
Human Swiss Prot No.:	O43521
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Bim (Phospho Thr56/116) Rabbit pAb at Human, Mouse, Rat
Storage and Stability:	Store at -20° C/1 year