

Antibodies Product Datasheet

Product Name:	Beclin-1 (Phospho Ser93/96) Antibody
Catalog Number:	C51140
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Beclin-1 (Phospho Ser93/96)
Concentration:	1 mg/ml
Gene Name:	BECN1 GT197
Protein Name:	Beclin-1 (Coiled-coil myosin-like BCL2-interacting protein) (Protein GT197)
Other Name:	Beclin-1 (Coiled-coil myosin-like BCL2-interacting protein) (Protein GT197)
Human Swiss Prot No.:	Q14457
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Beclin-1 (Phospho Ser93/96) Rabbit pAb at Human, Mouse, Rat
Storage and Stability:	Store at -20° C/1 year