

Antibodies Product Datasheet

Product Name:	GAB2 (Phospho Tyr614) Antibody
Catalog Number:	C51143
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GAB2 (Phospho Tyr614)
Concentration:	1 mg/ml
Gene Name:	GAB2 KIAA0571
Protein Name:	GRB2-associated-binding protein 2 (GRB2-associated binder 2) (Growth factor receptor bound protein 2-associated protein 2) (pp100)
Other Name:	GRB2-associated-binding protein 2 (GRB2-associated binder 2) (Growth factor receptor bound protein 2-associated protein 2) (pp100)
Human Swiss Prot No.:	Q9UQC2
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GAB2 (Phospho Tyr614) Rabbit pAb at Human, Mouse, Rat
Storage and Stability:	Store at -20° C/1 year