

Antibodies Product Datasheet

Product Name:	Androgen Receptor (Phospho Tyr267) Antibody
Catalog Number:	C51144
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Androgen Receptor (Phospho Tyr267)
Concentration:	1 mg/ml
Gene Name:	AR DHTR NR3C4
Protein Name:	Androgen receptor (Dihydrotestosterone receptor) (Nuclear receptor subfamily 3 group C member 4)
Other Name:	Androgen receptor (Dihydrotestosterone receptor) (Nuclear receptor subfamily 3 group C member 4)
Human Swiss Prot No.:	P10275
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Androgen Receptor (Phospho Tyr267) Rabbit pAb at Human, Mouse, Rat
Storage and Stability:	Store at -20° C/1 year