

Antibodies Product Datasheet

Product Name:	FGFR4 (Phospho Tyr754) Antibody
Catalog Number:	C51150
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Rabbit, polyclonal
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FGFR4 (Phospho Tyr754)
Concentration:	1 mg/ml
Gene Name:	FGFR4 JTK2 TKF
Protein Name:	Fibroblast growth factor receptor 4 (FGFR-4) (EC 2.7.10.1) (CD antigen CD334)
Other Name:	Fibroblast growth factor receptor 4 (FGFR-4) (EC 2.7.10.1) (CD antigen CD334)
Human Swiss Prot No.:	P22455
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FGFR4 (Phospho Tyr754) Rabbit pAb at Human, Mouse, Rat
Storage and Stability:	Store at -20° C/1 year