

Antibodies Product Datasheet

Product Name:	G6PC Antibody
Catalog Number:	C51166
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C-terminal G6PC
Concentration:	1 mg/ml
Gene Name:	G6PC G6PT
Protein Name:	Glucose-6-phosphatase (G-6-Pase) (G6Pase) (EC 3.1.3.9) (Glucose-6-phosphatase alpha) (G6Pase-alpha)
Other Name:	Glucose-6-phosphatase (G-6-Pase) (G6Pase) (EC 3.1.3.9) (Glucose-6-phosphatase alpha) (G6Pase-alpha)
Human Swiss Prot No.:	P35575
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of G6PC at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year