

Antibodies Product Datasheet

Product Name:	ITIH5 Antibody
Catalog Number:	C51182
Reactivity:	Human, Mouse
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C-terminal ITIH5
Concentration:	1 mg/ml
Gene Name:	ITIH5 KIAA1953 PP14776 UNQ311/PRO354
Protein Name:	Inter-alpha-trypsin inhibitor heavy chain H5 (ITI heavy chain H5) (ITI-HC5) (Inter-alpha-inhibitor heavy chain 5)
Other Name:	Inter-alpha-trypsin inhibitor heavy chain H5 (ITI heavy chain H5) (ITI-HC5) (Inter-alpha-inhibitor heavy chain 5)
Human Swiss Prot No.:	Q86UX2
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ITIH5 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year