

Antibodies Product Datasheet

Product Name:	CHCH1 Antibody
Catalog Number:	C51191
Reactivity:	Human, Mouse
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human N-terminal CHCH1
Concentration:	1 mg/ml
Gene Name:	CHCHD1 C10orf34
Protein Name:	Coiled-coil-helix-coiled-coil-helix domain-containing protein 1 (Nuclear protein C2360)
Other Name:	Coiled-coil-helix-coiled-coil-helix domain-containing protein 1 (Nuclear protein C2360)
Human Swiss Prot No.:	Q96BP2
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CHCH1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year