

Antibodies Product Datasheet

Product Name:	LIPT Antibody
Catalog Number:	C51209
Reactivity:	Human, Mouse
Applications:	IHC, WB
Dilutions:	WB 1:500-2000 IHC 1:50-200
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human N-terminal LIPT
Concentration:	1 mg/ml
Gene Name:	LIPT1
Protein Name:	Lipoyltransferase 1, mitochondrial (EC 2.3.1.-) (Lipoate biosynthesis protein) (Lipoate-protein ligase) (Lipoyl ligase)
Other Name:	Lipoyltransferase 1, mitochondrial (EC 2.3.1.-) (Lipoate biosynthesis protein) (Lipoate-protein ligase) (Lipoyl ligase)
Human Swiss Prot No.:	Q9Y234
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LIPT at Human, Mouse
Storage and Stability:	Store at -20° C/1 year